

# Gas Sampling Bags


## Gas Sampling Bags with PP Locking Combo Valve with Septum



Pull stem out to close



Push stem in to open



with stem  
pull out, turn  
stem 1/8"  
of a turn  
clockwise

Teflon Faced  
Septum

### Ordering Information

Description	QTY	Part No.
<b>ALTEF® Bags</b>		
6 x 6 inch, 0.6 L	1	061R-2101
7 x 7 inch, 1.0 L	1	061R-2102
9 x 9 inch, 2.0 L	1	061R-2103
10 x 10 inch, 3.0 L	1	061R-2104
12 x 12 inch, 5.0 L	1	061R-2105
12 x 17 inch, 8.0 L	1	061R-2106
12 x 18 inch, 9.5 L	1	061R-2107
12 x 19 inch, 10.0 L	1	061R-2108
12 x 21 inch, 12.0 L	1	061R-2109
18 x 18 inch, 16.0 L	1	061R-2110
18 x 24 inch, 25.0 L	1	061R-2111
24 x 24 inch, 40.0 L	1	061R-2112
24 x 30 inch, 56.0 L	1	061R-2113
24 x 36 inch, 73.0 L	1	061R-2114
30 x 30 inch, 80.0 L	1	061R-2115
30 x 36 inch, 100.0 L	1	061R-2116
<b>TEFLON® FEP Bags</b>		
6 x 6 inch, 0.6 L	1	061R-2201
7 x 7 inch, 1.0 L	1	061R-2202
6 x 10 inch, 1.2 L	1	061R-2203
6 x 12 inch, 1.5 L	1	061R-2204
9 x 9 inch, 2.0 L	1	061R-2205
10 x 10 inch, 3.0 L	1	061R-2206
12 x 12 inch, 5.0 L	1	061R-2207
12 x 18 inch, 9.5 L	1	061R-2208
12 x 19 inch, 10.0 L	1	061R-2209
18 x 18 inch, 16.0 L	1	061R-2210
18 x 24 inch, 25.0 L	1	061R-2211
24 x 24 inch, 40.0 L	1	061R-2212
24 x 30 inch, 56.0 L	1	061R-2213
24 x 36 inch, 73.0 L	1	061R-2214
30 x 30 inch, 80.0 L	1	061R-2215
<b>Multi-Layer Foil Bags</b>		
7 x 7 inch, 1.0 L	1	061R-2301
9 x 9 inch, 2.0 L	1	061R-2302
10 x 10 inch, 3.0 L	1	061R-2303
12 x 12 inch, 5.0 L	1	061R-2304
12 x 19 inch, 10.0 L	1	061R-2305
12 x 21 inch, 12.0 L	1	061R-2306