

## SGT Click-On Inline Traps

- Remove impurities to deliver 99.9999% pure gas
- No open gas line when changing the trap
- Diffusion-proof with glass and metal or Glass construction
- Pressure resistant up to 160psi and maximum flow 25L/min
- Estimated lifetime of more than 2 years

SGT Click-On Inline Traps are inline traps designed with Click-On adaptor connectors which allows inline cartridges to be exchanged without introducing oxygen. Spring loaded check valves seal when a filter is removed and open only when a new filter has been locked in place. Click-On Inline Traps measure 20 x 3.2cm.



Capacity				
Click-On Inline Traps	Moisture	Oxygen	Hydrocarbon <sup>4</sup>	Indicator/Color change
Moisture	21g	-	-	-
Hydrocarbon	-	-	36g	-
Oxygen	-	3000mL	-	-
Fuel <sup>1</sup>	10g	-	18g	-
Combi <sup>2</sup>	10g	1500mL	-	-
Triple <sup>3</sup>	6g	1000mL	12g	-
Triple <sup>3</sup> , He packed	6g	1000mL	12g	-
Triple <sup>3</sup> , H <sub>2</sub> packed	6g	1000mL	12g	-
Triple <sup>3</sup> , Indicator, He packed	0.1g	100mL	0.07g	H <sub>2</sub> O : Orange to clear, HC: Orange to green, O <sub>2</sub> : Green to gray

<sup>1</sup>Moisture/Hydrocarbon, <sup>2</sup>Moisture/Oxygen, <sup>3</sup>Moisture/Hydrocarbon/Oxygen, <sup>4</sup>As *n*-butane



Ordering Information		
Description	QTY	Part No.
<b>Click-On Inline Super-Clean™ Traps Kit with 1/8" Connector</b>		
Moisture Trap Kit	1	061P-3011
Oxygen Trap Kit	1	061P-3021
Hydrocarbon Trap Kit	1	061P-3031
Fuel Trap Kit	1	061P-3041
Combi Trap Kit	1	061P-3051
Triple Trap Kit	1	061P-3061
Triple, Helium packed Trap Kit	1	061P-3071
Triple, Hydrogen packed Trap Kit	1	061P-3081
Triple, Indicator, He packed Kit	1	061P-3091
<b>Replacement Traps</b>		
Moisture Trap	1	061P-3101
Oxygen Trap	1	061P-3111
Hydrocarbon Trap	1	061P-3121
Fuel Trap	1	061P-3131
Combi Trap	1	061P-3141
Triple Trap	1	061P-3151
Triple, Helium packed Trap	1	061P-3161
Triple, Hydrogen packed Trap	1	061P-3171
Triple, Indicator, He packed	1	061P-3181
<b>Connectors and Accessories</b>		
Brass Connectors, 1/8"	1	061P-3192
SS Connectors, 1/8"	1	061P-3202
SS Double Connector, 1/8"	1	061P-3211
Wall-Mounting Clamp	1	061P-3224
O-rings includes 10 large and 10 small	1	061P-3239