

Gas Purifiers

SGT Click-On Inline Indicating Traps

- Remove impurities to deliver 99.9999% pure gas
- No open gas line when changing the trap
- Diffusion-proof with glass and metal construction
- Pressure resistant up to 160psi and maximum flow 25L/min
- Estimated lifetime of more than 2 years

SGT Click-On Inline Indicating Traps are inline indicating traps designed with Click-On adaptor connectors which allows inline cartridges to be exchanged without introducing oxygen. Spring loaded check valves seal when a filter is removed and open only when a new filter has been locked in place. Click-On Inline Indicating Traps measure 20 x 3.2cm.



061P-4041



061P-4061

Capacity				
Click-On Inline Indicator Traps	Moisture	Oxygen	Hydrocarbon ⁴	Indicator/Color change
Moisture	6g	-	-	Orange to clear
Hydrocarbon	-	-	15g	Orange to green
Oxygen	-	850mL	-	Green to gray
Fuel ¹	3g	-	7g	H ₂ O : Orange to clear, HC: Orange to green
Combi ²	3g	415mL	-	H ₂ O : Orange to clear, O ₂ : Green to gray
Triple ³ , He packed	2g	400mL	1.5g	H ₂ O : Orange to clear, HC: Orange to green, O ₂ : Green to gray
Triple ³ , H ₂ packed ³	2g	400mL	1.5g	H ₂ O : Orange to clear, HC: Orange to green, O ₂ : Green to gray

¹Moisture/Hydrocarbon, ²Moisture/Oxygen, ³Moisture/Hydrocarbon/Oxygen, ⁴As *n*-butane

GC Accessories



061P-3192



061P-3211



061P-3224



061P-3239

Ordering Information		
Description	QTY	Part No.
Click-On Inline Indicator Trap Kits with 1/8" Connector		
Moisture Trap Kit	1	061P-4011
Oxygen Trap Kit	1	061P-4021
Hydrocarbon Trap Kit	1	061P-4031
Fuel Trap Kit	1	061P-4041
Combi Trap Kit	1	061P-4051
Triple, Helium packed Trap Kit	1	061P-4061
Triple, Hydrogen packed Trap	1	061P-4071
Replacement Traps		
Moisture Trap	1	061P-4081
Oxygen Trap	1	061P-4091
Hydrocarbon Trap	1	061P-4101
Fuel Trap	1	061P-4111
Combi Trap	1	061P-4121
Triple, Helium packed Trap	1	061P-4131
Triple, Hydrogen packed Trap	1	061P-4141
Connectors and Accessories		
Brass Connectors, 1/8"	2	061P-3192
SS Connectors, 1/8"	2	061P-3202
SS Double Connector, 1/8"	1	061P-3211
Wall-Mounting Clamp	4	061P-3224
O-rings includes 10 large and 10 small	20	061P-3239