

## Autosampler Syringes



064A-1011

064A-2031

### Ordering Information

Volume, Code, Gauge, Length, Point Style	QTY	OEM No.	PerkinElmer No.	Part No.
<b>Hamilton Syringes</b>				
5µL, ASN, 23s, 70mm, #3	1	88035	N6101390	24-88035
5µL, ASN, 26s, 70mm, #3	1	88040	N6101380	24-88040
<b>SGE Syringes</b>				
0.5µL, R, 23s, 70mm, Cone	1	000478	N6101252	21-000478
0.5µL, R, 26s, 70mm, Cone	1	000475		21-000475
5µL, F, 23s, 70mm, Cone	1	001954	N6101390	21-001954
5µL, F, 23s, GT, 70mm, Cone	1	001957		21-001957
5µL, F, 26s, 70mm, Cone	1	001953	N6101380	21-001953
5µL, F, 26s, GT, 70mm, Cone	1	001955		21-001955
50µL, F, 23s, 70mm, Cone	1	004670	N6101760	21-004670
Needles & Plunger 0.5µL, 23s, 70mm, Cone	1	037765	N6101469	21-037765
Needles & Plunger 0.5µL, 26s, 70mm, Cone	1	037750		21-037750

### Needle Termination Codes

#### Hamilton

ASN Autosampler  
Cemented Needle

#### SGE

F Fixed Needle  
R Removable Needle  
GT Gas-Tight



23s (0.64mm)



26s (0.47mm)



Point Style #3  
Blunt needle tip



Point Style Cone  
Cone tip