

## Autosampler Syringes




066A-8011

066A-9031

Needle Termination Codes	
<b>Hamilton</b>	
N	Cemented Needle
RN	Removable Needle
<b>SGE</b>	
F	Fixed Needle
R	Removable Needle
GT	Gas-Tight

Ordering Information				
Volume, Code, Gauge, Length, Point Style	QTY	OEM No.	Shimadzu No.	Part No.
<b>Hamilton Syringes for AOC-9</b>				
5µL, N, 26s, 51mm, #2	1	87900		24-87900
5µL, RN, 26s, 51mm, #2	1	87930		24-87930
10µL, RN, 26s, 51mm, #2	1	80330		24-80330
Needles, 26s, 51mm, #2	6	7758-02		24-7758-02
<b>SGE Syringes for AOC-14/17/20/20i</b>				
0.5µL, R, 23s, 42mm, Cone	1	000445		21-000445
0.5µL, R, 26s, 42mm, Cone	1	000440		21-000440
5µL, F, 23s, 42mm, Cone	1	001988		21-001988
5µL, F, 26s, 42mm, Cone	1	001987		21-001987
10µL, R, 23s, 42mm, Cone	1	002898		21-002898
10µL, R, 26s, 42mm, Cone	1	002897		21-002897
10µL, R, GT, 23s, 42mm, Cone	1	002902		21-002902
250µL, F, GT, 26s, 42mm, Cone	1	006682		21-006682
Needles, 23s, 42mm, Cone for 0.5µL syringes	1	033745		21-033745
Needles, 26s, 42mm, Cone for 0.5µL syringes	1	033738		21-033738
Needles, 23s, 42mm, Cone for 10µL syringes	2	037747		21-037747
Needles, 26s, 42mm, Cone for 10µL syringes	2	037745		21-037745



## Septum Nut



Septum nut provides consistent septum tightness, fewer bent needles, and easy leak detection. The nut contains a needle guide, ensuring that the needle consistently penetrates the septum in the same place, prolonging septum life. The guide also prevents the needle from striking the edge of the column or bending during insertion.

### Septum Nuts for Shimadzu GCs

- Ensure a leak-tight injection port, increase septum lifetime, and decrease maintenance requirements.
- Needle guide allows easy penetration and prevents septum coring.
- Manual injection septum nut allows use of 23-gauge needles for on-column injections.
- Made of high-quality stainless steel.

Ordering Information			
Description	QTY	Shimadzu No.	Part No.
<b>Septum Nut for Shimadzu 17A, 2010, and 2014 GCs</b>			
Septum Nut (for 23-Gauge Needles)	1	221-41286-00	066B-0011