

Liquid Phases

Ordering Information				
Description	Solvent*	Temp	QTY	Part No.
Liquid Phases				
Hi-EFF-2AP	C	100/210	25g	0900-0148
Hi-EFF-1BP	C	20/200	25g	0900-0149
Hi-EFF-2BP	C	100/200	25g	0900-0150
Igepal CO-630	M	100/200	50g	0900-0151
Igepal CO-880	C	170/200	50g	0900-0152
Igepal CO-990	C	100/220	50g	0900-0153
β,β,-Iminodipropionitrile	M	1100	25g	0900-0154
Kel-F Oil No. 10	A	20/100	50g	0900-0155
Lexan (polycarbonate resin)	C	220/270	25g	0900-0156
Mannitol		170/200	25g	0900-0157
Neopentyl Glycol Adipate	C	50/225	25g	0900-0158
Neopentyl Glycol Succinate	C	50/225	10g	0900-0159
Nujol (paraffin oil)	T	20/100	50g	0900-0160
β,β,-Oxydipropionitrile	M	20/100	25g	0900-0161
Phenyldiethanolamine	A	0/150	25g	0900-0162
Poly-A-103	C	70/275	10g	0900-0163
Poly-I-110	C	90/275	10g	0900-0164
Poly-S-176	C	150/400	3g	0900-0165
Poly-S-179	C	200/400	3g	0900-0166
Polyphenyl Ether (5 rings) OS-124	C	20/200	25g	0900-0167
Polypropylene Glycol	M	0/150	50g	0900-0168
Polyvinylpyrrolidinone (PVP)	M	20/200	50g	0900-0169
Propylene Glycol	C	0/50	50g	0900-0170
Reoplex 400 (polyester)	A	20/220	50g	0900-0171
Sebaconitrile	C	0/150	20g	0900-0172
Silicone DC-11	C	20/300	50g	0900-0173
Silicone DC-200, 350cstk	C	20/220	50g	0900-0174
Silicone DC-200, 12,500 cstk	C,T	0/200	50g	0900-0175
Silicone DC-410	C	20/250	50g	0900-0176
Silicone DC-401			50g	0900-0177
Silicone DC-550	C	20/225	50g	0900-0178
Silicone DC-704	A	20/250	50g	0900-0179
Silicone DC-710	A	20/225	50g	0900-0180
Silicone DC-HiVac Grease			50g	0900-0181
Silicone DC-QF-1	A	20/250	50g	0900-0182
Silicone DC-FS-1265	A	20/225	50g	0900-0183
Silicone GE-SE-30	C	50/300	50g	0900-0184
Silicone GE-SE-30 GC Grade	C	50/350	10g	0900-0185
Silicone GE-SE-52	C	50/300	50g	0900-0186
Silicone GE-SE-54	C	100/300	50g	0900-0187
Silicone GE-SF-96	C	20/250	50g	0900-0188
Silicone GE-XE-60	A	20/250	10g	0900-0189
Silicone GE-XF-1150		20/250	25g	0900-0190
Silicone OV-1	T	100/350	10g	0900-0191
Silicone OV-3	A	20/350	25g	0900-0192
Silicone OV-7	A	20/350	25g	0900-0193
Hi-EFF-1AP	A	0/350	20g	0900-0194

***Solvent Code:**

- A = Acetone
- C = Chloroform
- E = Ethyl Acetate
- T = Toluene
- M = Methanol