

# Liquid phases

Ordering Information				
Description	Solvent*	Temp	QTY	Part No.
<b>Liquid Phases</b>				
Silicone OV-17	A	20/350	25g	0900-0195
Silicone OV-22	A	20/350	10g	0900-0196
Silicone OV-25	A	20/350	10g	0900-0197
Silicone OV-61	A	20/350	10g	0900-0198
Silicone OV-73	T	20/350	10g	0900-0199
Silicone OV-101	T	20/350	20g	0900-0200
Silicone OV-105	A	120/250	10g	0900-0201
Silicone OV-202	C	20/275	10g	0900-0202
Silicone OV-210	C	20/275	25g	0900-0203
Silicone OV-215	E	20/275	10g	0900-0204
Silicone OV-225	A	20/275	10g	0900-0205
Silicone OV-275	A	20/275	5g	0900-0206
Silicone OV-330	A	30/250	5g	0900-0207
OV-351 (Replaces FFAP)	C	50/250	10g	0900-0208
Silicone OV-1701	A	20/325	3g	0900-0209
Silicone Silar 5CP	C	50/275	5g	0900-0210
Silicone Silar 7CP	C	50/275	5g	0900-0211
Silicone Silar 9CP	C	50/275	5g	0900-0212
Silicone Silar 10C	C	50/275	5g	0900-0213
SP-2100	T	20/350	50g	0900-0214
SP-2250	A	20/350	50g	0900-0215
SP-2401	C	0/275	50g	0900-0216
SP-1000	C	50/250	50g	0900-0217
SP-2300	C,A	50/275	50g	0900-0218
SP-2310	A	50/275	50g	0900-0219
SP-2330	A	50/275	50g	0900-0220
SP-2340	A	50/275	50g	0900-0221
Silicone UCC L-45			50g	0900-0222
Silicone UCC W-98		0/300	50g	0900-0223
Span 80 (sorbitan monooleate)	T	20/150	25g	0900-0224
Squalene	T	20/150	50g	0900-0225
Sucrose Acetate Isobutyrate (SAIB)	T	30/200	25g	0900-0226
Tergitol NPX	C	10/175	50g	0900-0227
THEED	M	20/125	25g	0900-0228
Tricresyl Phosphate	M	20/125	50g	0900-0229
Triethanolamine	M	25/75	50g	0900-0230
1,2,3-Tris(2-cyanoethoxy)propane (TCEP)	M	29/150	50g	0900-0231
Triton X-100	A	20/190	50g	0900-0232
Triton X-305	A	20/250	50g	0900-0233
Tween 80	M	20/160	50g	0900-0234
UCON LB-550-X	M	20/200	50g	0900-0235
UCON 50-HB-280-X	M	20/200	50g	0900-0236
UCON 50-HB-2000	A	20/200	50g	0900-0237
UCON 50-HB-5100	M	20/200	50g	0900-0238
Versamid 900	A	190/275	50g	0900-0239

**\*Solvent Code:**

- A = Acetone
- C = Chloroform
- E = Ethyl Acetate
- T = Toluene
- M = Methanol