


Screw Thread Headspace Vials

Screw Thread Headspace Vials

- Available in clear, Type I Class A or amber, Type I Class B Borosilicate glass with 18mm screw thread
- Not require crimper nor decapper
- Multiple turn threading offer a tight seal
- Choice of 10 mL and 20mL
- Round bottom vials are compatible with most autosampler
- Pre-assembled magnetic caps and septa minimize contamination
- Compatible with Waters®, Shimadzu, Perkin Elmer, Jasco and Varian autosampler



Ordering Information

Description	QTY	Part No.
18mm Screw Thread Headspace Vials		
10mL Clear Glass, Round Bottom	100	0709-1011
10mL Amber Glass, Round Bottom	100	0709-1021
20mL Clear Glass, Round Bottom	100	0709-1031
20mL Amber Glass, Round Bottom	100	0709-1041
18mm Magnetic Caps and Septa		
Cap with 0.115" Natural PTFE/Blue Silicone	100	0709-4011
Cap with 0.045" White PTFE/Blue Silicone	100	0709-4021
Cap with 0.120" White PTFE/Chlorobutyl	100	0709-4031
Septa		
Blue PTFE/White Silicone	100	0709-3011