


Membrane Filters and Accessories

Glass Syringes

- Reinforced flanges with two flat sides to prevent rolling
- All borosilicate glass plunger and barrel for accurate, leak-proof performance and resistance to chemicals and thermal shock.
- Graduations are fused on to the syringe barrel
- Interchangeable plungers



5mL Glass syringe with glass luer slip tip



20mL Glass syringe with metal luer lock tip

Plastic Syringes

- Made of medical-grade polypropylene
- Sterile Disposable
- Ideal for use with syringe filters and general purpose applications
- Come in an easy-to-open sealed package
- Eliminate Ion Contamination in IC Applications



Filtration

Ordering Information		
Size	QTY	Part No.
Glass Syringe, Luer-Tip		
2mL	1	0211-1211
	12	0211-1213
5mL	1	0211-1411
	12	0211-1413
10mL	1	0211-1511
	12	0211-1513
20mL	1	0211-1611
	12	0211-1613
50mL	1	0211-1811
	12	0211-1813
Glass Syringe, Luer-Lock		
2mL	1	0211-2211
	12	0211-2213
5mL	1	0211-2411
	12	0211-2413
10mL	1	0211-2511
	12	0211-2513
20mL	1	0211-2611
	12	0211-2613
50mL	1	0211-2811
	12	0211-2813

Ordering Information		
Size	QTY	Part No.
Plastic Syringe, Luer-Tip		
1mL	100	0211-1125
3mL	100	0211-1325
5mL	100	0211-1425
10mL	100	0211-1525
20mL	50	0211-1624
Plastic Syringe, Luer-Lock		
3mL	100	0211-2325
5mL	100	0211-2425
10mL	100	0211-2525
20mL	50	0211-2624