



# Fittings

## SealTight™ Two-Piece Fittings Two-Piece Fittings

- Upchurch Scientific® SealTight™ Two-Piece Fingertight, pressure rated to 9,000 psi
- All SealTight Nuts are for use with 1/16" OD tubing, 10-32 port and are designed to be used exclusively with the SealTight™ Ferrule
- Upchurch Scientific®'s original Two-Piece Fingertight
- Works with 1/16" OD Tubing and female 10-32 ports
- Fingertight without tools
- Pressure rated to 6,000 psi



HPLC Accessories

## SS Two-Piece Fittings

- 316 Stainless Steel Fittings are rated to 10,000 psi when wrench tightened
- Works with 1/16" OD Tubing and female 10-32 ports



### Ordering Information

Description	QTY	Part No.
<b>Two-Piece Fittings with Ferrules</b>		
Derlin Wing Nuts, Red	10	05AC-2017
PEEK Wing Nuts, Natural	10	05AC-2027
PEEK Long Nuts, Natural	10	05AC-2037
PEEK Short Nuts, Natural	10	05AC-2047
PEEK Ferrules, Natural	10	05AC-2057
<b>SealTight™ Two-Piece Fittings with Ferrules</b>		
PEEK Short Nuts, Hex Head, Black	10	05AC-2067
PEEK Long Nuts, Hex Head, Black	10	05AC-2077
PEEK Short Nuts, Standard Head, Black	10	05AC-2087
PEEK Long Nuts, Standard Head, Black	10	05AC-2097
PEEK Short Nuts, Headless, Black	10	05AC-2107
PEEK Long Nuts, Headless, Black	10	05AC-2117
PEEK Long Nuts, Knurl/Hex Head, Black	10	05AC-2127
PEEK SealTight™ Ferrules, Black	10	05AC-2137
<b>SS Fittings with Ferrules</b>		
Standard Nuts	10	05AC-2147
Standard Ferrules	10	05AC-2157
Waters® Nuts	10	05AC-2167
Waters® Ferrules	10	05AC-2177
SSI Nuts	10	05AC-2187
SSI Ferrules	10	05AC-2297
Rheodyne® Nuts with Ferrules	10	05AC-2207
Rheodyne® Long Nuts with Ferrules	10	05AC-2217
Rheodyne® Extra-Long Nuts with Ferrules	10	05AC-2227
Rheodyne® Ferrules	10	05AC-2237