

## SGE Syringes for Autosamplers



05BG-1011

05GG-1031



05IG-1011

05KG-1011

Ordering Information			
Description	QTY	OEM No.	Part No.
<b>Agilent 1090, 1100</b>			
SGE 25D-HP1090, 25µL GASTIGHT, 1/4-32 UNEF	1	003670	21-003670
SGE 250D-HP1090, 250µL GASTIGHT, 1/4-32 UNEF	1	006670	21-006670
Replacement plunger for Syringe 25µL	1	0318181	21-0318181
Replacement plunger for Syringe 250µL	1	031829	21-031829
<b>Perkin-Elmer L-200</b>			
SGE 50D-CX, 50µL GASTIGHT, 1/4-28 UNF	1	004995	21-004995
SGE 100D-CX-GT 100µL GASTIGHT, 1/4-28 UNF	1	005990	21-005990
SGE 250D-CX-GT 250µL GASTIGHT, 1/4-28 UNF	1	006995	21-006995
SGE 500D-CX-GT 500µL GASTIGHT, 1/4-28 UNF	1	007995	21-007995
SGE 1MD-C-GT 1mL GASTIGHT, 1/4-28 UNF	1	008185	21-008185
SGE 2.5MD-C-GT 2.5mL GASTIGHT, 1/4-28 UNF	1	008687	21-008687
Replacement plunger for 50µL	1	0318221	21-0318221
Replacement plunger for 100µL	1	0318271	21-0318271
Replacement plunger for 250µL	1	0318341	21-0318341
Replacement plunger for 500µL	1	0318381	21-0318381
Replacement plunger for 1mL	1	0318441	21-0318441
Replacement plunger for 2.5mL	1	031854	21-031854
<b>Thermo Finnigan, AS300, AS 100, AS3000, AS1000</b>			
SGE 250C-TS 250µL gas tight, 1/4-28 UNF	1	006660	21-006660
SGE 500C-TS 500µL gas tight, 1/4-28 UNF	1	007680	21-007680
SGE 1MC-TS 1mL gas tight, 1/4-28 UNF	1	008180	21-008180
SGE 2.5MTS 2.5mL gas tight, 1/4-28 UNF	1	008660	21-008660
SGE 5MC-TS 5mL gas tight, 1/4-28 UNF	1	008780	21-008780
Replacement plunger for 250µL	1	031833	21-031833
Replacement plunger for 500µL	1	031838	21-031838
Replacement plunger for 1mL	1	031844	21-031844
Replacement plunger for 2.5mL	1	031853	21-031853
Replacement plunger for 5mL	1	031856	21-031856
<b>Waters/WISP/Alliance</b>			
SGE 25D-WISP 25µL, 1/4-28 UNF	1	003990	21-003990
SGE 250D-WISP 250µL, 1/4-28 UNF	1	006690	21-006690
Replacement plunger for 25µL	1	031819	21-031819
Replacement plunger for 250µL	1	031834	21-031834
WISP Dipper Needle	1	038265	21-038265

## Syringe Cleaners



Ordering Information		
Description	QTY	Part No.
<b>Syringe Cleaners</b>		
120 VAC	1	24-76610
220 VAC	1	24-76615
Replacement Septa, 5mm Diameter	12	24-76620A
<b>Needle Cleaning Kit</b>		
Kit, Includes: cleaning wires (listed below) and 70mL of Cleaning concentrate, 500mL	1	24-76620
<b>Cleaning Wires, 7" Length</b>		
.00300" OD, for Gauge 23s, 26s, 31-33	1	24-18306
.00350" OD for Gauge 26s, 31-33	10	24-18300
.00497" OD for Gauge 22s, 25s, 28-30	10	24-18301
.00659" OD for Gauge 27	10	24-18302
.00815" OD for Gauge 24-26	10	24-18303
.01207" OD for Gauge 22, 23	10	24-18304