

HPLC Syringes

Hamilton Syringes for Manual Injection



Syringes for Rheodyne® Injection Valves

Ordering Information					
Volume, Termination ¹ , Gauge, Length, PS#	QTY	OEM Model	OEM No.	Part No.	
Hamilton Fixed Needle Syringes					
10µL, N, 22s, 51mm, 3	1	701SNR	80365	24-80365	
25µL, N, 22s, 51mm, 3	1	702SNR	80465	24-80465	
50µL, N, 22s, 51mm, 3	1	705SNR	80565	24-80565	
100µL, N, 22s, 51mm, 3	1	710SNR	80665	24-80665	
250µL, N, 22s, 51mm, 3	1	725SNR	80765	24-80765	
500µL, N, 22s, 51mm, 3	1	750SNR	80865	24-80865	
Hamilton Fixed Needle, Small Bore, GASTIGHT Syringes					
25µL, N, 22s, 51mm, 3	1	1702N	80275	24-80275	
50µL, N, 22s, 51mm, 3	1	1705N	80975	24-80975	
100µL, N, 22s, 51mm, 3	1	1710N	81075	24-81075	
250µL, N, 22s, 51mm, 3	1	1725N	81175	24-81175	
Hamilton Fixed Needle, Large Bore, GASTIGHT Syringes					
25µL, N, 22, 51mm, 3	1	1702N	80285	24-80285	
50µL, N, 22, 51mm, 3	1	1705N	80985	24-80985	
100µL, N, 22, 51mm, 3	1	1710N	81085	24-81085	
250µL, N, 22, 51mm, 3	1	1725N	81185	24-81185	
500µL, N, 22, 51mm, 3	1	1750LTN	81216	24-81216	
1mL, N, 22, 51mm, 3	1	1001LTN	81316	24-81316	
2.5mL, N, 22, 51mm, 3	1	1002LTN	81416	24-81416	
5.0mL, N, 22, 51mm, 3	1	1005LTN	81516	24-81516	
Hamilton Removable Needle, GASTIGHT Syringes					
10µL, RN, 22s, 51mm, 3	1	1701RNR	80065	24-80065	
25µL, RN, 22s, 51mm, 3	1	1702RNR	80265	24-80265	
50µL, RN, 22s, 51mm, 3	1	1705RNR	80965	24-80965	
100µL, RN, 22s, 51mm, 3	1	1710RNR	81065	24-81065	
250µL, RN, 22s, 51mm, 3	1	1725RNR	81165	24-81165	
500µL, RN, 22s, 51mm, 3	1	1750RSNR	81265	24-81265	
1mL, RN, 22s, 51mm, 3	1	1001RNR	81365	24-81365	
Replacement Needles					
10-100µL RN Needles, 22s, 51mm, 3	6		7770-01	24-7770-01	
10-100µL RN Needles, 22, 51mm, 3	6		7770-02	24-7770-02	
250µL-5.0mL RN Needles, 22s, 51mm, 3	6		7780-03	24-7780-03	
250µL-5.0mL RN Needles, 22, 51mm, 3	6		7780-04	24-7780-04	



Point Style #3

Blunt needle tip

Blunt needle point for use with HPLC injection valves and for sample pipetting

Needle Gauge

Gauge	Normal OD (mm)	Normal ID (mm)
25s	0.51	0.25
25s	0.72	0.15
22	0.72	0.41

Syringes for Waters® U6K Injection Valves

Ordering Information					
Volume, Termination ¹ , Gauge, Length, PS#	QTY	OEM Model	OEM No.	Part No.	
Hamilton Removable Needle, GASTIGHT Syringes					
25µL, N, 25s, 50mm, 3	1	1702RNW	80238	24-80238	
50µL, N, 25s, 50mm, 3	1	1705RNW	80938	24-80938	
100µL, N, 25s, 50mm, 3	1	1710RNW	81038	24-81038	
250µL, N, 25s, 50mm, 3	1	1725RNW	81138	24-81138	
Hamilton Removable Needle, GASTIGHT Syringes with Reinforce Plunger					
25µL, N, 25s, 50mm, 3	1	1802RNW	84980	24-84980	
50µL, N, 25s, 50mm, 3	1	1805RNW	84983	24-84983	
100µL, N, 25s, 50mm, 3	1	1810RNW	84986	24-84986	
250µL, N, 25s, 50mm, 3	1	1825RNW	84989	24-84989	
Replacement Needles					
RN Needles, 25s, 50mm, 3	1		8647-01	24-8647-01	