


# Tubings

## Stainless Steel Tubings

- Presize cut 316 SS tubing
- Square , burr-free ends



Ordering Information		
Description	QTY	Part No.
<b>0.005" x 1/16" (Red)</b>		
5cm	1	05AD-2011
10cm	1	05AD-2021
20cm	1	05AD-2031
30cm	1	05AD-2041
50cm	1	05AD-2051
1m	1	05AD-2061
5feet	1	05AD-2071
25feet	1	05AD-2081
<b>0.007" x 1/16" (Black)</b>		
5cm	1	05AD-2091
10cm	1	05AD-2101
20cm	1	05AD-2111
30cm	1	05AD-2121
50cm	1	05AD-2131
1m	1	05AD-2141
5feet	1	05AD-2151
25feet	1	05AD-2161
<b>0.010" x 1/16" (Blue)</b>		
5cm	1	05AD-2171
10cm	1	05AD-2181
20cm	1	05AD-2191
30cm	1	05AD-2201
50cm	1	05AD-2211
1m	1	05AD-2221
5feet	1	05AD-2231
25feet	1	05AD-2241
<b>0.020" x 1/16" (Yellow)</b>		
5cm	1	05AD-2251
10cm	1	05AD-2261
20cm	1	05AD-2271
30cm	1	05AD-2281
50cm	1	05AD-2291
1m	1	05AD-2301
5feet	1	05AD-2311
25feet	1	05AD-2321

Ordering Information		
Description	QTY	Part No.
<b>0.030" x 1/16" (White)</b>		
5cm	1	05AD-2331
10cm	1	05AD-2341
20cm	1	05AD-2351
30cm	1	05AD-2361
50cm	1	05AD-2371
1m	1	05AD-2381
5feet	1	05AD-2391
25feet	1	05AD-2401
<b>0.040" x 1/16"</b>		
5cm	1	05AD-2411
10cm	1	05AD-2421
20cm	1	05AD-2431
30cm	1	05AD-2441
50cm	1	05AD-2451
1m	1	05AD-2461
5feet	1	05AD-2471
25feet	1	05AD-2481
<b>0.046" x 1/16"</b>		
5cm	1	05AD-2491
10cm	1	05AD-2501
20cm	1	05AD-2511
30cm	1	05AD-2521
50cm	1	05AD-2531
1m	1	05AD-2541
5feet	1	05AD-2551
25feet	1	05AD-2561
<b>0.080" x 1/16"</b>		
15cm	1	05AD-2571
25cm	1	05AD-2581
1m	1	05AD-2591
3m	1	05AD-2601
5m	1	05AD-2611

## Stainless Steel Tubings

- These SS tubing kits include 10 pieces – 5cm long,
- 10 pieces – 10cm long

Ordering Information		
Description	QTY	Part No.
<b>SS Tubing Kits</b>		
Tubing kit 0.010" ID	1	05AD-2621
Tubing kit 0.007" ID	1	05AD-2631
Tubing kit 0.005" ID	1	05AD-2641
Tubing kit 0.020" ID	1	05AD-2651
Assorted 1/16" OD tubing kit, includes 10 x 05AD-2251, 10 x 05AD-2261, 5 x 05AD-2271, 10 x 05AD-2171, 10 x 05AD-2181 and 5 x 05AD-2191	1	05AD-2661