

HPLC Instrument Parts

Thermo: LDC/Milton Roy HPLC Systems



05IH-0041



05IH-0051



05IH-0151



05IH-0071



05IH-0111



05IH-0061

Ordering Information			
Description	QTY	Thermo No.	Part No.
Replacement parts for Thermo: LDC/Milton Roy HPLC Systems			
Inlet check valve assembly	1	900947001	05IH-0011
Outlet check valve assembly	1	900947002	05IH-0021
Replacement check valve cartridge	1		05IH-0032
1/8" Nut with Ferrule for Suction Tube	1	1960066976	05IH-0041
1/16" Nut with Ferrule for Discharge Tube	1	1960066975	05IH-0051
Sapphire piston	1	801306	05IH-0061
Seal, UHMW-PE	1	206234	05IH-0071
Seal UHMW-PE Replacement Kit	1	801598	05IH-0081
Seal, UHMW-PE, Gold	1	206234	05IH-0091
Seal UHMW-PE Gold Replacement Kit	1	801892001	05IH-0101
Seals, Black Teflon	1	206129001	05IH-0111
Seal Black Teflon Replacement Kit	1	801598	05IH-0121
Crosshead Spring, Constromatic II/Mini-pump	1	916208000	05IH-0131
Crosshead Spring, Minimetric	1	204040	05IH-0141
D2 lamp for 3100/3200/4000/4100/4200	1		05IH-0151
D2 lamp for Thermo UV 6000 (DAD)	1	108052	05IH-0301

- Black Teflon seals are PTFE with a carbon fill, this seal is a very good general purpose seal for any solvent.
- Gold plunger seals are made of UHMW-PE (Ultra High Molecular Weight Polyethylene) with fill that gives it a gold appearance. The seal provides good seal life and does not shed particles like that of other seals. The gold seal is inert to all common HPLC solvent mixes.

Thermo: Spectra-Physics HPLC Systems



05IH-0161



05IH-0171



05IH-0281



05IH-0211



05IH-0201

Ordering Information			
Description	QTY	Thermo No.	Part No.
Replacement parts for Thermo: Spectra-Physics HPLC			
Inlet check valve assmb. 8700/8800/8810/ISOCHROM/P	1	A3495-010	05IH-0161
Outlet check valve assmb. 8700/8800/8810/ISOCHROM/P	1	A3490-010	05IH-0171
Transducer valve assmb. 8700/8800/8810/ISOCHROM/P	1	A3499-010	05IH-0181
Replacement check valve cartridge, 2/Pk	1		05IH-0191
Sapphire Piston	1	A3102-010	05IH-0201
Seal, UHMW-PE for 8800/8810/ISOCHROM/P	1	A2962-010	05IH-0211
Seal, Black Teflon for 8800/8810/ISOCHROM/P	1	A2962-010	05IH-0221
Back-up Seal, Kel-F for 8800/8810/ISOCHROM/P	1	A2973-010	05IH-0231
Push Seal, UHMW-PE for 8800/8810/ISOCHROM/P	1	A2963-010	05IH-0241
Push Seal, Black Teflon for 8800/8810/ISOCHROM/P	1	A2963-010	05IH-0251
Seal, UHMW-PE for 8700/8700XR/8750	1	A1703-020	05IH-0261
Seal, Black Teflon for 8700/8700XR/8750	1	A1703-020	05IH-0271
D2 lamp for Focus	1	9551-0023	05IH-0281
D2 Lamp Hight Intensity	1	9551-0202	05IH-0291

- Black Teflon seals are PTFE with a carbon fill, this seal is a very good general purpose seal for any solvent.
- Gold plunger seals are made of UHMW-PE (Ultra High Molecular Weight Polyethylene) with fill that gives it a gold appearance. The seal provides good seal life and does not shed particles like that of other seals. The gold seal is inert to all common HPLC solvent mixes.