

HPLC Instrument Parts

Varian, Gilson/Rainin HPLC Systems



05JH-0121



05JH-0131



05JH-0141

Ordering Information			
Description	QTY	Varian No.	Part No.
Replacement parts for Varian, Gilson/Rainin HPLC Systems			
Inlet check valve, Sapphire for Rainin 5,10&25ML Heads	1	R007101616	05JH-0011
Inlet check valve, Ceramic for Rainin 5,10&25ML Heads	1	R007101678	05JH-0011
Outlet check valve, Sapphire for Rainin 5,10&25ML Heads	1	R007101617	05JH-0021
Outlet check valve, Ceramic for Rainin 5,10&25ML Heads	1	R007101679	05JH-0021
Replacement check valve cartridge	2		05JH-0032
Inlet check valve, Sapphire for Varian Prostar 210/215	1	R007101616	05JH-0011
Inlet check valve, Ceramic for Varian Prostar 210/215	1	R007101678	05JH-0011
Outlet check valve, Sapphire for Varian Prostar 210/215	1	R007101617	05JH-0021
Outlet check valve, Ceramic for Varian Prostar 210/215	1	R007101679	05JH-0021
Replacement check valve cartridge	2		05JH-0032
Inlet check valve assembly for Varian 2010	1	00-997261-09	05JH-0041
Outlet check valve assembly for Varian 2010	1	00-997261-10	05JH-0051
Replacement check valve cartridge	2		05JH-0062
Inlet check valve assembly for Varian 2510	1	00-997554-19	05JH-0071
Outlet check valve assembly for Varian 2510	1	00-997554-18	05JH-0081
Replacement check valve cartridge	2		05JH-0092
Outlet check valve assembly for Varian 9010/9012	1	03-919465-90	05JH-0101
Replacement check valve cartridge	2		05JH-0112
Seal, UHMW-PE for Gilson 10mL S/SC heads	1	5463125895	05JH-0121
Seal, UHMW-PE for Rainin 10mL S/SC heads	1	R007101637	05JH-0121
Seal, Black Teflon for Gilson 10mL S/SC heads	1	5463125095	05JH-0131
Seal, Black Teflon for Rainin 10mL S/SC heads	1	R007101636	05JH-0131
Piston for Gilson 10mL S/SC heads	1	365009	05JH-0141
Piston for Rainin 10mL S/SC heads	1	R007101658	05JH-0141
Seal, UHMW-PE for Gilson 5mL S/SC heads	1	5463088863	05JH-0151
Seal, UHMW-PE for Rainin 5mL S/SC heads	1	R007101634	05JH-0151
Seal, Black Teflon for Gilson 5mL S/SC heads	1	5463088863	05JH-0161
Seal, Black Teflon for Rainin 5mL S/SC heads	1	R007101633	05JH-0161
Piston for Gilson 5mL S/SC heads	1	365008	05JH-0171
Piston for Rainin 5mL S/SC heads	1	R007101657	05JH-0181
Seal, UHMW-PE for Varian 2010/2510	1	00-997261-37	05JH-0191
Seal, Black Teflon for Varian 2010/2510	1	00-997261-37	05JH-0201
Piston for Varian 2010/2510	1	00-997261-08	05JH-0211
Seal, UHMW-PE for Varian 9012/9012	1	03-919397-00	05JH-0221
Seal, Black Teflon for Varian 9012/9012	1	03-919397-00	05JH-0231
Piston for Varian 2010/2510	1	03-919102-90	05JH-0241
Hydraulic Module with Damper, Purge Valve and Outlet Manifold	1	212-LC	05JH-0311
D2 lamp for Gilson/Rainin 112	1		05JH-0251
D2 Lamp for Gilson/Rainin DAD	1		05JH-0261
D2 lamp for Varian 2550	1		05JH-0271
D2 Lamp for Varian 2050	1		05JH-0281
D2 lamp for Varian Prostar 330 DAD	1		05JH-0291
D2 Lamp for Varian Prostar 310/UV-50/100/200/9050	1		05JH-0301

- Black Teflon seals are PTFE with a carbon fill, this seal is a very good general purpose seal for any solvent.
- Gold plunger seals are made of UHMW-PE (Ultra High Molecular Weight Polyethylene) with fill that gives it a gold appearance. The seal provides good seal life and does not shed particles like that of other seals. The gold seal is inert to all common HPLC solvent mixes.