


12x32mm Crimp Top Wide Mouth Vials

12x32mm Crimp Top Wide Mouth Vials

- Available in clear, Type I Class A or amber, Type I Class B Borosilicate glass or high quality polypropylene
- Total Volume of 2mL and useable volume of 1.5mL
- Glass vials are available with or without spot
- Wide mouth with Crimp Top neck
- Aluminum or Polypropylene open-top crimp seals with choices of septa
- Limited volume inserts are available for small volume samples
- A vial kit in PP Gridded Box includes vials, Caps with septa
- Compatible with Agilent, Waters®, HTA, Shimadzu, Thermo and Varian autosampler



Ordering Information

Description	QTY	Part No.
12x32mm Crimp Top Wide Mouth Vials		
Clear glass	100	0704-1011
Clear glass, spot	100	0704-1021
Amber glass	100	0704-1041
Amber glass, spot	100	0704-1051
Crimp Seals with Septa		
Al seals, PTFE/Red Rubber	100	0703-4011
Al seals, PTFE/White Silicone	100	0703-4021
Al seals, PTFE/Silicone/PTFE	100	0703-4031
Al seals, PTFE/Silicone, Slit	100	0703-4041
Vial Kits in Gridded Box		
Clear glass spot vials with Al seals and PTFE/Red Rubber Septa	100	0704-5011
Amber glass spot vials with Al seals and PTFE/Red Rubber Septa	100	0704-5021
Clear glass spot vials with Al seals and PTFE/White Silicone Septa	100	0704-5031
Amber glass spot vials with Al seals and PTFE/White Silicone Septa	100	0704-5041
Inserts		
Glass, 100µL, w/Bottom Spring	100	0701-6011
Glass, 250µL, w/Bottom Spring	100	0701-6021
Glass, 300µL	100	0701-6031
Glass, 400µL, Flat Bottom	100	0701-6041