

## EPA Vials



- Available in clear, Type I Class A or amber, Type I Class B Borosilicate glass with 24/400 screw thread finish
- For EPA CFR Part 136 an EPA 40 CFR 141 for Volatile organic analysis
- Choice of 20, 30, 40 and 60mL
- Vial kit includes vials, caps and 0.125" PTFE/White Silicone septa or top-hat septa
- Open-Top caps are white or black polypropylene for sampling through the septa
- Solid-Top caps are white polypropylene for sample storage

### Ordering Information

Description	QTY	Part No.
<b>EPA Vial Kit</b>		
20mL <sup>1</sup> Clear, Black cap and septa	72	0712-1011
20mL Amber, Black cap and septa	72	0712-1021
20mL Clear, Black cap and Top-hat septa	72	0712-1031
20mL Amber, Black cap and Top-hat septa	72	0712-1041
20mL Clear, White Solid-Top cap and septa	72	0712-1051
20mL Amber, White Solid-Top cap and septa	72	0712-1061
30mL <sup>2</sup> Clear, White cap and septa	100	0712-1071
30mL Amber, White cap and septa	100	0712-1081
30mL Clear, Black cap and Top-hat septa	100	0712-1091
30mL Amber, Black cap and Top-hat septa	100	0712-1101
30mL Clear, White Solid-Top cap and septa	100	0712-1111
30mL Amber, White Solid-Top cap and septa	100	0712-1121
40mL <sup>3</sup> Clear, White cap and septa	100	0712-1131
40mL Amber, White cap and septa	100	0712-1141
40mL Clear, Black cap and Top-hat septa	72	0712-1151
40mL Amber, Black cap and Top-hat septa	72	0712-1161
40mL Clear, White Solid-Top cap and septa	100	0712-1171
40mL Amber, White Solid-Top cap and septa	100	0712-1181
60mL <sup>4</sup> Clear, White cap and septa	100	0712-1191
60mL Amber, White cap and septa	100	0712-1201
60mL Clear, Black cap and Top-hat septa	100	0712-1211
60mL Amber, Black cap and Top-hat septa	100	0712-1221
60mL Clear, White Solid-Top cap and septa	100	0712-1231
60mL Amber, White Solid-Top cap and septa	100	0712-1241
<b>Caps</b>		
Black Open-Top PP Cap	100	0712-2011
White Solid-Top PP Cap with PTFE Septa	100	0712-2021
White Solid-Top PP Cap with PTFE/Silicone Septa	100	0712-2031
<b>Septa</b>		
0.060" PTFE/Silicone Septa	100	0712-3011
0.125" PTFE/Silicone Septa	100	0712-3021
<b>Caps with Septa</b>		
White Open-Top PP Cap, 0.060" PTFE/Silicone	100	0712-4011
White Open-Top PP Cap, 0.125" PTFE/Silicone	100	0712-4021

<sup>1</sup> 20mL vial measure 28x57mm with 24/400 Finish

<sup>2</sup> 30mL vial measure 28x73mm with 24/400 Finish

<sup>3</sup> 40mL vial measure 28x95mm with 24/400 Finish

<sup>4</sup> 60mL vial measure 28x140mm with 24/400 Finish