


8x40mm Shell Vials

8x40mm Shell Vials

- Available in clear, Type I Class A or amber, Type I Class B Borosilicate glass or high quality polypropylene
- Usable Volume of 1mL for Flat Bottom style and 700µL for PP limited volume
- Polyethylene Snap plug with Starburst for easy needle penetration
- Limited volume inserts are available for small volume samples
- A kit includes shell vials and polyethylene snap plug with Starburst
- For Waters® WISP™ 96 Style autosampler



| Ordering Information | | |
|--|-----|-----------|
| Description | QTY | Part No. |
| 8x40mm Shell Vials | | |
| Clear glass, Flat Bottom, 1mL | 100 | 0710-1011 |
| Amber glass, Flat Bottom, 1mL | 100 | 0710-1021 |
| Snap plug with Starburst, 8mm | | |
| Clear PE Plug | 100 | 0710-2011 |
| 8x40mm Shell Vials kit with plugs | | |
| Clear glass, Flat Bottom, 1mL | 100 | 0710-5011 |
| Amber glass, Flat Bottom, 1mL | 100 | 0710-5021 |
| Polypropylene, Flat Bottom, 1mL | 250 | 0710-5032 |
| Polypropylene, Limited Volume, 700µL | 250 | 0710-5042 |
| Inserts | | |
| Glass, 250µL, w/Bottom Spring | 100 | 0701-6021 |
| Glass, 300µL | 100 | 0701-6031 |
| Polypropylene, 300µL | 100 | 0701-6051 |